


Tropomi SWIR ISRF [laser 1]

coverage : 2020-01-03T03:24:01
duration : 303
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T01:20:07


normalized response (flattened)



Tropomi SWIR ISRF [laser 1]

coverage : 2020-01-03T03:24:01
duration : 303
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T01:20:07

normalized response (flattened)



Tropomi SWIR ISRF [laser 1]

coverage : 2020-01-03T03:24:01
duration : 303
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T01:20:07

normalized response

