


Tropomi SWIR ISRF [laser 1]

coverage : 2020-04-13T15:30:33
duration : 304
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T02:31:05


normalized response (flattened)



Tropomi SWIR ISRF [laser 1]

coverage : 2020-04-13T15:30:33
duration : 304
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T02:31:05

normalized response (flattened)



Tropomi SWIR ISRF [laser 1]

coverage : 2020-04-13T15:30:33
duration : 304
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T02:31:06

normalized response

