


Tropomi SWIR ISRF [laser 1]

coverage : 2020-07-24T03:29:45
duration : 303
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T03:42:13


normalized response (flattened)



Tropomi SWIR ISRF [laser 1]

coverage : 2020-07-24T03:29:45
duration : 303
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T03:42:13

normalized response (flattened)



Tropomi SWIR ISRF [laser 1]

coverage : 2020-07-24T03:29:45
duration : 303
column : [470 471]
peltier current : [10.0, 12.5]%
diffuser : no oscillation
created : 2023-08-24T03:42:13

normalized response

